


# Sling Sizing Guide

## Sizes

Use this sizing chart as a general guide. Keep in mind that two residents/patients that are the same weight may have a completely different body shape and size and may require different size slings.

Weight	Sling Size	Stand Belt Size	Strap Color Code
>100 lbs.	Small	Small	Black
100-210 lbs.	Medium	Medium	Blue
210-310 lbs.	Large	Large	Yellow
310-450 lbs.	X-Large	X-Large	White
450-1000 lbs.	XX-Large		Purple



## Medicare Slings Contain Antimicrobial Product Protection

Toileting, incontinence, wound care and simple everyday use make slings and belts very susceptible to microbe growth. It is also common practice for slings and belts to be used multiple times on multiple people during the day before they are washed.

Microban protection provides continuous antimicrobial cleaning that inhibits the growth of microbes such as bacteria, mold and mildew that can cause stains, odors and product deterioration. Built-in during the manufacturing process, Microban antimicrobial protection will not wash or wear off, providing protection for the useful lifetime of the product and keeping the slings and belts cleaner between cleanings.

