

Sling & Belt Sizing Guide



Sling Sizing

Use the sizing chart below as a general guide. Keep in mind that two patients that are the same weight may have a completely different body shape and size and may require different size slings.

Weight	Sling Size	Strap Color
0–100 lbs.	Small	Black
100–210 lbs.	Medium	Blue
210–310 lbs.	Large	Yellow
310–450 lbs.	X-Large	White
450–1000 lbs.	Bariatric	Purple



Belt Sizing

Use this sizing chart as a general guide. Keep in mind that two patients that are the same weight may have a completely different body shape and size and may require different size belts.

Belt Size	Belt Waist Size	Strap Color
Small	24" to 48"	Black
Medium	30" to 54"	Blue
Large	36" to 60"	Yellow
X-Large	42" to 66"	White



Medcare Slings & Belts Contain Antimicrobial Protection

Toileting, incontinence, wound care and simple everyday use make slings and belts very susceptible to microbe growth. It is also common practice for slings and belts to be used multiple times on multiple people during the day before they are washed.

Antimicrobial protection works 24/7 for an added level of cleanliness protection. Antimicrobial protection inhibits the growth of microbes such as bacteria, mold and mildew that can cause stains, odors and product deterioration.

Built-in during the manufacturing process, antimicrobial protection will not wash or wear off, providing protection for the useful lifetime of the product and keeping the slings and belts cleaner between cleanings.